

# Recombinant Human CMC4

Catalog No: #GP11769



Package Size: #GP11769-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CMC4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	C-x(9)-C motif containing 4
Other Names	p8; C6.1B; MTCP1; MTCP1B; MTCP1NB; p8MTCP1
Accession No.	Swissprot:P56277 Gene Accession:BC002600
Uniprot	P56277
GeneID	100272147;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene was identified by involvement in some t(X;14) translocations associated with mature T-cell proliferations. This region has a complex gene structure, with a common promoter and 5' exon spliced to two different sets of 3' exons that encode two different proteins. This gene represents the downstream 8 kDa protein that localizes to mitochondria.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only