

Recombinant Human MTFR1

Catalog No: #GP11771



Package Size: #GP11771-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MTFR1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 179-304 amino acids of human mitochondrial fission regulator 1
Target Name	mitochondrial fission regulator 1
Other Names	CHPPR; FAM54A2
Accession No.	Swissprot:Q15390Gene Accession:BC043498
Uniprot	Q15390
GeneID	9650;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a mitochondrial protein that is characterized by a poly-proline rich region. A chicken homolog of this protein promotes mitochondrial fission and the mouse homolog protects cells from oxidative stress. A related pseudogene of this gene is found on chromosome X.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only