

Recombinant Human NDUFA9

Catalog No: #GP11786



Package Size: #GP11786-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NDUFA9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 36-300 amino acids of human NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 9, 39kDa
Target Name	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 9, 39kDa
Other Names	CC6; CI39k; CI-39k; NDUFS2L; SDR22E1
Accession No.	Swissprot:Q16795Gene Accession:BC015837
Uniprot	Q16795
GenID	4704;
Storage	-20~-80°C, pH 7.6 PBS

Background

The encoded protein is a subunit of the hydrophobic protein fraction of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. A pseudogene has been identified on chromosome 12.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only