

## Recombinant Human PVRL2

Catalog No: #GP11803



Package Size: #GP11803-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PVRL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 32-300 amino acids of human poliovirus receptor-related 2 (herpesvirus entry mediator B)
Target Name	poliovirus receptor-related 2 (herpesvirus entry mediator B)
Other Names	HVEB; PRR2; CD112; PVRR2
Accession No.	Swissprot:Q92692Gene Accession:BC003091
Uniprot	Q92692
GenID	5819;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only