

# Recombinant Human PVRL3

Catalog No: #GP11805



Package Size: #GP11805-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Recombinant Human PVRL3   |
| Brief Description     | Recombinant Protein   |
| Immunogen Description | Fusion protein corresponding to C terminal 124 amino acids of human poliovirus receptor-related 3 |
| Target Name           | poliovirus receptor-related 3   |
| Other Names           | PPR3; PRR3; CD113; PVRR3; CDW113; NECTIN-3  |
| Accession No.         | Swissprot:Q9NQS3Gene Accession:BC067808   |
| Uniprot               | Q9NQS3  |
| GeneID                | 25945;  |
| Storage               | -20~-80°C, pH 7.6 PBS   |

## Background

This gene encodes a member of the nectin family of proteins, which function as adhesion molecules at adherens junctions. This family member interacts with other nectin-like proteins and with afadin, a filamentous actin-binding protein involved in the regulation of directional motility, cell proliferation and survival. This gene plays a role in ocular development involving the ciliary body. Mutations in this gene are believed to result in congenital ocular defects. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only