

Recombinant Human NLGN1

Catalog No: #GP11807



Package Size: #GP11807-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human NLGN1 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 44-300 amino acids of human neuroligin 1 |
| Target Name | neuroligin 1 |
| Other Names | NL1 |
| Accession No. | Swissprot:Q8N2Q7Gene Accession:BC032555 |
| Uniprot | Q8N2Q7 |
| GeneID | 22871; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only