

Recombinant Human NCSTN

Catalog No: #GP11808



Package Size: #GP11808-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NCSTN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 34-330 amino acids of human nicastrin
Target Name	nicastrin
Other Names	ATAG1874
Accession No.	Swissprot:Q92542Gene Accession:BC047621
Uniprot	Q92542
GeneID	23385;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a type I transmembrane glycoprotein that is an integral component of the multimeric gamma-secretase complex. The encoded protein cleaves integral membrane proteins, including Notch receptors and beta-amyloid precursor protein, and may be a stabilizing cofactor required for gamma-secretase complex assembly.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only