

Recombinant Human PITPNM1

Catalog No: #GP11809



Package Size: #GP11809-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PITPNM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human phosphatidylinositol transfer protein, membrane-associated 1
Target Name	phosphatidylinositol transfer protein, membrane-associated 1
Other Names	Rd9; NIR2; RDGB; DRES9; RDGB1; RDGBA; PITPNM; RDGBA1
Accession No.	Swissprot:O00562Gene Accession:BC022230
Uniprot	O00562
GeneID	9600;
Storage	-20~-80°C, pH 7.6 PBS

Background

PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the *Drosophila* retinal degeneration B (rdgB) protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only