

Recombinant Human NIT1

Catalog No: #GP11810



Package Size: #GP11810-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NIT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 47-243 amino acids of human nitrilase 1
Target Name	nitrilase 1
Accession No.	Swissprot:Q86X76Gene Accession:BC046149
Uniprot	Q86X76
GeneID	4817;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the nitrilase protein family with homology to bacterial and plant nitrilases, enzymes that cleave nitriles and organic amides to the corresponding carboxylic acids plus ammonia. Multiple transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only