

Recombinant Human NPM3

Catalog No: #GP11817



Package Size: #GP11817-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NPM3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 11-145 amino acids of human nucleophosmin/nucleoplasmin 3
Target Name	nucleophosmin/nucleoplasmin 3
Other Names	PORMIN; TMEM123
Accession No.	Swissprot:O75607Gene Accession:BC054868
Uniprot	O75607
GenID	10360;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is related to the nuclear chaperone phosphoproteins, nucleoplasmin and nucleophosmin. This protein is strongly expressed in diverse cell types where it localizes primarily to the nucleus. Based on its similarity to nucleoplasmin and nucleophosmin, this protein likely functions as a molecular chaperone in the cell nucleus.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only