

Recombinant Human HMGN5

Catalog No: #GP11820



Package Size: #GP11820-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HMGN5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 130 amino acids of human high mobility group nucleosome binding domain 5
Target Name	high mobility group nucleosome binding domain 5
Other Names	NSBP1; NBP-45
Accession No.	Swissprot:P82970Gene Accession:BC005342
Uniprot	P82970
GenID	79366;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a nuclear protein with similarities to the high mobility group proteins, HMG14 and HMG17, which suggests that this protein may function as a nucleosomal binding and transcriptional activating protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only