

## Recombinant Human NT5C1A

Catalog No: #GP11821



Package Size: #GP11821-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NT5C1A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human 5'-nucleotidase, cytosolic IA
Target Name	5'-nucleotidase, cytosolic IA
Other Names	CN1; CNI; CN-I; CN1A; CN-IA
Accession No.	Swissprot:Q9BXI3 Gene Accession:BC103880
Uniprot	Q9BXI3
GeneID	84618;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Cytosolic nucleotidases, such as NT5C1A, dephosphorylate nucleoside monophosphates. Dephosphorylates the 5' and 2'(3')-phosphates of deoxyribonucleotides and has a broad substrate specificity. Helps to regulate adenosine levels in heart during ischemia and hypoxia.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only