

Recombinant Human GNL3

Catalog No: #GP11823



Package Size: #GP11823-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GNL3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 282-456 amino acids of human guanine nucleotide binding protein-like 3
Target Name	guanine nucleotide binding protein-like 3
Other Names	NS; E2IG3; NNP47; C77032
Accession No.	Swissprot:Q9BVP2Gene Accession:BC001024
Uniprot	Q9BVP2
GenelD	26354;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene may interact with p53 and may be involved in tumorigenesis. The encoded protein also appears to be important for stem cell proliferation. This protein is found in both the nucleus and nucleolus. Three transcript variants encoding two different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only