

Recombinant Human PAGE5

Catalog No: #GP11841



Package Size: #GP11841-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAGE5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	P antigen family, member 5 (prostate associated)
Other Names	CT16; CT16.1; CT16.2; GAGEE1; PAGE-5
Accession No.	Swissprot:Q96GU1Gene Accession:BC009230
Uniprot	Q96GU1
GeneID	90737;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the GAGE family, which is expressed in a variety of tumors and in some fetal and reproductive tissues. The protein encoded by this gene shares a sequence similarity with other GAGE/PAGE proteins. It may also belong to a family of CT (cancer-testis) antigens. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene, but the biological validity of some variants have not been determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only