

Recombinant Human PAGE1

Catalog No: #GP11842



Package Size: #GP11842-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PAGE1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	P antigen family, member 1 (prostate associated)
Other Names	AL5; CT16.3; GAGE-9; GAGEB1; PAGE-1
Accession No.	Swissprot:O75459Gene Accession:BC004861
Uniprot	O75459
GeneID	8712;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to a family of genes that are expressed in a variety of tumors but not in normal tissues, except for the testis. Unlike the other gene family members, this gene does not encode an antigenic peptide. Nothing is presently known about the function of this protein.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only