

# Recombinant Human SCGB2A1

Catalog No: #GP11850



Package Size: #GP11850-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SCGB2A1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-95 amino acids of human SCGB2A1
Target Name	secretoglobin, family 2A, member 1
Other Names	LPHC; LPNC; MGB2; UGB3
Accession No.	Swissprot:O75556Gene Accession:BC062218
Uniprot	O75556
GeneID	4246;
Storage	-20~-80°C, pH 7.6 PBS

## Background

SCGB2A1 may bind androgens and other steroids, may also bind estramustine, a chemotherapeutic agent used for prostate cancer. It may be under transcriptional regulation of steroid hormones.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only