

## Recombinant Human PDGFRL

Catalog No: #GP11852



Package Size: #GP11852-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PDGFRL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 272-375 amino acids of human platelet-derived growth factor receptor-like
Target Name	platelet-derived growth factor receptor-like
Other Names	PDGRL; PRLTS
Accession No.	Swissprot:Q15198Gene Accession:BC010927
Uniprot	Q15198
GenID	5157;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein with significant sequence similarity to the ligand binding domain of platelet-derived growth factor receptor beta. Mutations in this gene, or deletion of a chromosomal segment containing this gene, are associated with sporadic hepatocellular carcinomas, colorectal cancers, and non-small cell lung cancers. This suggests this gene product may function as a tumor suppressor.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only