

Recombinant Human PDLIM1

Catalog No: #GP11853



Package Size: #GP11853-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDLIM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 257 amino acids of human PDZ and LIM domain 1
Target Name	PDZ and LIM domain 1
Other Names	CLIM1; CLP36; CLP-36; hCLIM1; HEL-S-112
Accession No.	Swissprot:O00151Gene Accession:BC000915
Uniprot	O00151
GeneID	9124;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the enigma protein family. The protein contains two protein interacting domains, a PDZ domain at the amino terminal end and one to three LIM domains at the carboxyl terminal. It is a cytoplasmic protein associated with the cytoskeleton. The protein may function as an adapter to bring other LIM-interacting proteins to the cytoskeleton. Pseudogenes associated with this gene are located on chromosomes 3, 14 and 17.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only