

Recombinant Human PEF1

Catalog No: #GP11857



Package Size: #GP11857-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PEF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 114-284 amino acids of human penta-EF-hand domain containing 1
Target Name	penta-EF-hand domain containing 1
Other Names	ABP32; PEF1A
Accession No.	Swissprot:Q9UBV8Gene Accession:BC012561
Uniprot	Q9UBV8
GenID	553115;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. The encoded protein has been shown to form a heterodimer with the programmed cell death 6 gene product and may modulate its function in Ca(2+) signaling. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 1.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only