

Recombinant Human PGRMC1

Catalog No: #GP11860



Package Size: #GP11860-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PGRMC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 44-195 amino acids of human PGRMC1
Target Name	progesterone receptor membrane component 1
Other Names	MPR; HPR6.6
Accession No.	Swissprot:O00264Gene Accession:BC034238
Uniprot	O00264
GeneID	10857;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a putative membrane-associated progesterone steroid receptor. The protein is expressed predominantly in the liver and kidney.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only