

Recombinant Human ANP32A

Catalog No: #GP11862



Package Size: #GP11862-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ANP32A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	acidic nuclear phosphoprotein 32 family member A
Other Names	LANP; MAPM; PP32; HPPCn; PHAP1; PHAPI; I1PP2A; C15orf1
Accession No.	Swissprot:P39687Gene Accession:BC125143
Uniprot	P39687
GeneID	8125;
Storage	-20~-80°C, pH 7.6 PBS

Background

Implicated in a number of cellular processes, including proliferation, differentiation, caspase-dependent and caspase-independent apoptosis, suppression of transformation (tumor suppressor), inhibition of protein phosphatase 2A, regulation of mRNA trafficking and stability in association with ELAVL1, and inhibition of acetyltransferases as part of the INHAT (inhibitor of histone acetyltransferases) complex. Plays a role in E4F1-mediated transcriptional repression.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only