

Recombinant Human PLK4

Catalog No: #GP11872



Package Size: #GP11872-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PLK4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 12-265 amino acids of human PLK4
Target Name	polo-like kinase 4
Other Names	SAK; STK18; MCCR2
Accession No.	Swissprot:O00444Gene Accession:BC036023
Uniprot	O00444
GeneID	10733;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the polo family of serine/threonine protein kinases. The protein localizes to centrioles, complex microtubule-based structures found in centrosomes, and regulates centriole duplication during the cell cycle. Three alternatively spliced transcript variants that encode different protein isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only