

Recombinant Human PNLIPRP1

Catalog No: #GP11873



Package Size: #GP11873-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PNLIPRP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 218-467 amino acids of human PNLIPRP1
Target Name	pancreatic lipase-related protein 1
Other Names	PLRP1
Accession No.	Swissprot:P54315Gene Accession:BC025784
Uniprot	P54315
GeneID	5407;
Storage	-20~-80°C, pH 7.6 PBS

Background

May function as inhibitor of dietary triglyceride digestion. Lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters (in vitro) (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only