

Recombinant Human PLTP

Catalog No: #GP11874



Package Size: #GP11874-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PLTP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human phospholipid transfer protein
Target Name	phospholipid transfer protein
Other Names	BPIFE; HDLCQ9
Accession No.	Swissprot:P55058Gene Accession:BC019847
Uniprot	P55058
GeneID	5360;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only