

Recombinant Human PMP2

Catalog No: #GP11876



Package Size: #GP11876-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PMP2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	peripheral myelin protein 2
Other Names	P2; MP2; FABP8; M-FABP
Accession No.	Swissprot:P02689Gene Accession:BC034997
Uniprot	P02689
GeneID	5375;
Storage	-20~-80°C, pH 7.6 PBS

Background

Peripheral Myelin protein-2 (PMP2, Myelin P2, MP2), a small basic enzyme, is one of the principal proteins of peripheral Myelin and presumably participates in the transport of fatty acids or the metabolism of Myelin lipids. Myelin P2 is similar in amino acid sequence and tertiary structure to fatty acid binding proteins found in the liver, adipocytes and intestine, although its expression is restricted to the nervous system.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only