

Recombinant Human PPP2R3C

Catalog No: #GP11886



Package Size: #GP11886-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PPP2R3C
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 180 amino acids of human protein phosphatase 2, regulatory subunit B", gamma
Target Name	protein phosphatase 2, regulatory subunit B", gamma
Other Names	G4-1; G5pr; C14orf10
Accession No.	Swissprot:Q969Q6Gene Accession:BC006823
Uniprot	Q969Q6
GeneID	55012;
Storage	-20~-80°C, pH 7.6 PBS

Background

Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit gamma?is an enzyme?that in humans is encoded by the?PPP2R3C?gene. May regulate MCM3AP phosphorylation through phosphatase recruitment. May play a role in the activation-induced cell death of B-cells By similarity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only