

Recombinant Human RMDN3

Catalog No: #GP11897



Package Size: #GP11897-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RMDN3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human regulator of microtubule dynamics 3
Target Name	regulator of microtubule dynamics 3
Other Names	RMD3; RMD-3; FAM82C; FAM82A2; ptpip51
Accession No.	Swissprot:Q96TC7Gene Accession:BC063844
Uniprot	Q96TC7
GeneID	55177;
Storage	-20~-80°C, pH 7.6 PBS

Background

Involved in cellular calcium homeostasis regulation. May participate in differentiation and apoptosis of keratinocytes. Overexpression induces apoptosis. Present at high level in epidermis and seminiferous epithelium: while basal cells in the epidermis and spermatogonia show no perceptible amount, keratinocytes of suprabasal layers and differentiating first-order spermatocytes up to spermatids exhibit high expression.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only