

Recombinant Human PTPN3

Catalog No: #GP11902



Package Size: #GP11902-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PTPN3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 29-312 amino acids of human protein tyrosine phosphatase, non-receptor type 3
Target Name	protein tyrosine phosphatase, non-receptor type 3
Other Names	PTPH1; PTP-H1
Accession No.	Swissprot:P26045Gene Accession:BC126117
Uniprot	P26045
GenID	5774;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This protein contains a C-terminal PTP domain and an N-terminal domain homologous to the band 4.1 superfamily of cytoskeletal-associated proteins.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only