

# Recombinant Human RAB25

Catalog No: #GP11909



Package Size: #GP11909-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RAB25
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 210 amino acids of human RAB25, member RAS oncogene family
Target Name	RAB25, member RAS oncogene family
Other Names	CATX-8; RAB11C
Accession No.	Swissprot:P57735Gene Accession:BC009831
Uniprot	P57735
GeneID	57111;
Storage	-20~-80°C, pH 7.6 PBS

## Background

RAB proteins, such as RAB25, are members of the RAS superfamily of small GTPases that are involved in membrane trafficking. Members of the RAB11 subfamily, including RAB25, control the return of internalized membrane-associated moieties to the cell surface.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only