

# Recombinant Human RAB2A

Catalog No: #GP11910



Package Size: #GP11910-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RAB2A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB2A, member RAS oncogene family
Other Names	LHX; RAB2
Accession No.	Swissprot:P61019Gene Accession:BC008929
Uniprot	P61019
GeneID	5862;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene belongs to the Rab family, members of which are small molecular weight guanosine triphosphatases (GTPases) that contain highly conserved domains involved in GTP binding and hydrolysis. The Rabs are membrane-bound proteins, involved in vesicular fusion and trafficking. This protein is a resident of pre-Golgi intermediates, and is required for protein transport from the endoplasmic reticulum (ER) to the Golgi complex.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only