

Recombinant Human RAB27A

Catalog No: #GP11911



Package Size: #GP11911-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB27A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB27A, member RAS oncogene family
Other Names	GS2; RAM; RAB27; HsT18676
Accession No.	Swissprot:P51159Gene Accession:BC132800
Uniprot	P51159
GeneID	5873;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only