

Recombinant Human RAB4A

Catalog No: #GP11912



Package Size: #GP11912-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB4A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB4A, member RAS oncogene family
Other Names	RAB4; HRES1; HRES-1; HRES-1/RAB4
Accession No.	Swissprot:P20338Gene Accession:BC002438
Uniprot	P20338
GeneID	5867;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only