

Recombinant Human RAB20

Catalog No: #GP11915



Package Size: #GP11915-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB20
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-234 amino acids of human RAB20
Target Name	RAB20, member RAS oncogene family
Accession No.	Swissprot:Q9NX57Gene Accession:BC026025
Uniprot	Q9NX57
GeneID	55647;
Storage	-20~-80°C, pH 7.6 PBS

Background

Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as *S.aureus* and *M.tuberculosis*. Plays a role in the fusion of phagosomes with lysosomes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only