

Recombinant Human RASAL1

Catalog No: #GP11921



Package Size: #GP11921-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RASAL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 300 amino acids of human RAS protein activator like 1 (GAP1 like)
Target Name	RAS protein activator like 1 (GAP1 like)
Other Names	RASAL
Accession No.	Swissprot:O95294Gene Accession:BC014420
Uniprot	O95294
GeneID	8437;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is member of the GAP1 family of GTPase-activating proteins. These proteins stimulate the GTPase activity of normal RAS p21 but not its oncogenic counterpart. Acting as a suppressor of RAS function, the protein enhances the weak intrinsic GTPase activity of RAS proteins resulting in the inactive GDP-bound form of RAS, thereby allowing control of cellular proliferation and differentiation.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only