

# Recombinant Human RASAL2

Catalog No: #GP11922



Package Size: #GP11922-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RASAL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human RAS protein activator like 2
Target Name	RAS protein activator like 2
Other Names	NGAP
Accession No.	Swissprot:Q9UJF2Gene Accession:BC110611
Uniprot	Q9UJF2
GeneID	9462;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein that contains the GAP-related domain (GRD), a characteristic domain of GTPase-activating proteins (GAPs). GAPs function as activators of Ras superfamily of small GTPases. The protein encoded by this gene is able to complement the defective RasGAP function in a yeast system. Two alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only