

# Recombinant Human RASSF4

Catalog No: #GP11924



Package Size: #GP11924-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RASSF4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 127-321 amino acids of human Ras association (RalGDS/AF-6) domain family member 4
Target Name	Ras association (RalGDS/AF-6) domain family member 4
Other Names	AD037
Accession No.	Swissprot:Q9H2L5?Gene Accession:BC032593
Uniprot	Q9H2L5
GeneID	83937;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only