

Recombinant Human RASSF8

Catalog No: #GP11925



Package Size: #GP11925-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RASSF8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 181 amino acids of human RASSF8
Target Name	Ras association (RalGDS/AF-6) domain family (N-terminal) member 8
Other Names	HOJ1; C12orf2
Accession No.	Swissprot:Q8NHQ8Gene Accession:BC030021
Uniprot	Q8NHQ8
GeneID	11228;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the Ras-association domain family (RASSF) of tumor suppressor proteins. This gene is essential for maintaining adherens junction function in epithelial cells and has a role in epithelial cell migration. It is a lung tumor suppressor gene candidate. A chromosomal translocation t(12;22)(p11.2;q13.3) leading to the fusion of this gene and the FBLN1 gene is found in a complex type of synpolydactyly. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only