

Recombinant Human RASSF9

Catalog No: #GP11926



Package Size: #GP11926-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RASSF9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human Ras association (RalGDS/AF-6) domain family (N-terminal) member 9
Target Name	Ras association (RalGDS/AF-6) domain family (N-terminal) member 9
Other Names	PAMCI; PCIP1; P-CIP1
Accession No.	Swissprot:O75901Gene Accession:BC031589
Uniprot	O75901
GeneID	9182;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene localizes to perinuclear endosomes. This protein associates with peptidylglycine alpha-amidating monooxygenase, and may be involved with the trafficking of this enzyme through secretory or endosomal pathways.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only