

## Recombinant Human RBBP5

Catalog No: #GP11931



Package Size: #GP11931-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RBBP5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 332-538 amino acids of human retinoblastoma binding protein 5
Target Name	retinoblastoma binding protein 5
Other Names	RBQ3; SWD1
Accession No.	Swissprot:Q15291Gene Accession:BC053856
Uniprot	Q15291
GenID	5929;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only