

# Recombinant Human RDH11

Catalog No: #GP11934



Package Size: #GP11934-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RDH11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 252 amino acids of human retinol dehydrogenase 11 (all-trans/9-cis/11-cis)
Target Name	retinol dehydrogenase 11 (all-trans/9-cis/11-cis)
Other Names	MDT1; CGI82; PSDR1; RALR1; SCALD; ARSDR1; HCBP12; SDR7C1
Accession No.	Swissprot:Q8TC12Gene Accession:BC011727
Uniprot	Q8TC12
GeneID	51109;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is an NADPH-dependent retinal reductase and a short-chain dehydrogenase/reductase. The encoded protein has no steroid dehydrogenase activity.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only