

Recombinant Human REG4

Catalog No: #GP11937



Package Size: #GP11937-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human REG4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 102 amino acids of human regenerating islet-derived family, member 4
Target Name	regenerating islet-derived family, member 4
Other Names	GISP; RELP; REG-IV
Accession No.	Swissprot:Q9BYZ8Gene Accession:BC017089
Uniprot	Q9BYZ8
GeneID	83998;
Storage	-20~-80°C, pH 7.6 PBS

Background

Calcium-independent lectin displaying mannose-binding specificity and able to maintain carbohydrate recognition activity in an acidic environment. May be involved in inflammatory and metaplastic responses of the gastrointestinal epithelium.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only