

Recombinant Human RERG

Catalog No: #GP11940



Package Size: #GP11940-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human RERG |
| Brief Description | Recombinant Protein |
| Immunogen Description | Full length fusion protein |
| Target Name | RAS-like, estrogen-regulated, growth inhibitor |
| Accession No. | Swissprot:Q96A58Gene Accession:BC007997 |
| Uniprot | Q96A58 |
| GeneID | 85004; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

RERG, a member of the RAS superfamily of GTPases, inhibits cell proliferation and tumor formation (Finlin et al., 2001 [PubMed 11533059]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only