

Recombinant Human RNASET2

Catalog No: #GP11950



Package Size: #GP11950-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RNASET2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 232 amino acids of human ribonuclease T2
Target Name	ribonuclease T2
Other Names	RNASE6PL; bA514O12.3
Accession No.	Swissprot:O00584Gene Accession:BC001660
Uniprot	O00584
GeneID	8635;
Storage	-20~-80°C, pH 7.6 PBS

Background

This ribonuclease gene is a novel member of the Rh/T2/S-glycoprotein class of extracellular ribonucleases. It is a single copy gene that maps to 6q27, a region associated with human malignancies and chromosomal rearrangement.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only