

Recombinant Human PDLIM4

Catalog No: #GP11952



Package Size: #GP11952-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDLIM4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	PDZ and LIM domain 4
Other Names	RIL
Accession No.	Swissprot:P50479Gene Accession:BC016765
Uniprot	P50479
GeneID	8572;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only