

Recombinant Human DSTYK

Catalog No: #GP11954



Package Size: #GP11954-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DSTYK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 635-884 amino acids of human DSTYK
Target Name	dual serine/threonine and tyrosine protein kinase
Other Names	RIP5; RIPK5; CAKUT1; DustyPK; HDCMD38P
Accession No.	Swissprot:Q6XUX3Gene Accession:BC053627
Uniprot	Q6XUX3
GeneID	25778;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a dual serine/threonine and tyrosine protein kinase which is expressed in multiple tissues. It is thought to function as a regulator of cell death. Multiple transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only