

## Recombinant Human RNF11

Catalog No: #GP11959

Package Size: #GP11959-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RNF11
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ring finger protein 11
Other Names	CGI-123; SID1669
Accession No.	Swissprot:Q9Y3C5Gene Accession:BC020964
Uniprot	Q9Y3C5
GeneID	26994;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene contains a RING-H2 finger motif, which is known to be important for protein-protein interactions. The expression of this gene has been shown to be induced by mutant RET proteins (MEN2A/MEN2B). The germline mutations in RET gene are known to be responsible for the development of multiple endocrine neoplasia (MEN).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only