

Recombinant Human RRAGC

Catalog No: #GP11974



Package Size: #GP11974-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RRAGC
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Ras-related GTP binding C
Other Names	GTR2; RAGC; TIB929
Accession No.	Swissprot:Q9HB90Gene Accession:BC016668
Uniprot	Q9HB90
GeneID	64121;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the GTR/RAG GTP-binding protein family. The encoded protein is a monomeric guanine nucleotide-binding protein which forms a heterodimer with RRAGA and RRAGB and is primarily localized to the cytoplasm. The encoded protein promotes intracellular localization of the mTOR complex. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only