

Recombinant Human MAPK12

Catalog No: #GP11989



Package Size: #GP11989-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MAPK12
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	mitogen-activated protein kinase 12
Other Names	ERK3; ERK6; SAPK3; PRKM12; SAPK-3; P38GAMMA
Accession No.	Swissprot:P53778Gene Accession:BC015741
Uniprot	P53778
GeneID	6300;
Storage	-20~-80°C, pH 7.6 PBS

Background

Activation of members of the mitogen-activated protein kinase family is a major mechanism for transduction of extracellular signals. Stress-activated protein kinases are one subclass of MAP kinases. The protein encoded by this gene functions as a signal transducer during differentiation of myoblasts to myotubes.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only