

Recombinant Human SDF4

Catalog No: #GP11993



Package Size: #GP11993-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SDF4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 312 amino acids of human stromal cell derived factor 4
Target Name	stromal cell derived factor 4
Other Names	Cab45; SDF-4
Accession No.	Swissprot:Q9BRK5Gene Accession:BC006211
Uniprot	Q9BRK5
GeneID	51150;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a stromal cell derived factor that is a member of the CREC protein family. The encoded protein contains six EF-hand motifs and calcium-binding motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only