

Recombinant Human SERTAD1

Catalog No: #GP11995



Package Size: #GP11995-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SERTAD1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	SERTA domain containing 1
Other Names	SEI1; TRIP-Br1
Accession No.	Swissprot:Q9UHV2Gene Accession:BC002670
Uniprot	Q9UHV2
GeneID	29950;
Storage	-20~-80°C, pH 7.6 PBS

Background

Acts at E2F-responsive promoters as coregulator to integrate signals provided by PHD- and/or bromodomain-containing transcription factors. Stimulates E2F1/TFDP1 transcriptional activity. Renders the activity of cyclin D1/CDK4 resistant to the inhibitory effects of CDKN2A/p16INK4A.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only